Références/Referencias/ Referencias/ Referências:	Composition du coffret/ Kit composition/ Composición del kit/ Conteúdo da embalagem :
ALSL-0250 8 x 25 mL	R1 8 x 20 mL + R2 8 x 5 mL
ALSL-0455 4 x 55 mL	R1 4 x 44 mL + R2 4 x 11 mL
ALSL-0410 2 x 62,5 mL	R1 2 x 50 mL + R2 1 x 26 mL
ALSL-0430 4 x 62,5 mL	R1 4 x 50 mL + R2 2 x 26 mL
ALSL-0510 5 x 125 mL	R1 5 x 100 mL + R2 1 x 127 mL

Références/Referencias/ Referencias/ Referências:	Composition du coffret/ Kit composition/ Composición del kit/ Conteúdo da embalagem :
ALSL-0250 8 x 25 mL	R1 8 x 20 mL + R2 8 x 5 mL
ALSL-0455 4 x 55 mL	R1 4 x 44 mL + R2 4 x 11 mL
ALSL-0410 2 x 62,5 mL	R1 2 x 50 mL + R2 1 x 26 mL
ALSL-0430 4 x 62,5 mL	R1 4 x 50 mL + R2 2 x 26 mL
ALSL-0510 5 x 125 mL	R1 5 x 100 mL + R2 1 x 127 mL

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BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Panteghini, M., Bais, R., *Enzyme, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
- Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders Company), (2006), 64
- Dufour, R., *The liver: Function and chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 586, appendix
- Schumann, G., et al. *Clin Chem Lab Med.*, (2002), **40**, 718.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations*, WHO/DIL/LAB/99.1 Rev.2, (2002), 21.
- Ward, M.K., Cockayne, S., *Enzymology. Clinical Chemistry : Concepts and Application*, Anderson, S.C., Cockayne, S. (W.B.Saunders eds. Philadelphia USA), (1993), 238.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
- Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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English - EN

INTENDED USE

ELITech Clinical Systems ALT/GPT 4+1 SL is intended for the quantitative *in vitro* diagnostic determination of Alanine amino transferase (ALT) in human serum and plasma.

CLINICAL SIGNIFICANCE (1-3,6)

Alanine aminotransferase (ALT), also known as glutamate pyruvate transaminase (GPT), is a transaminase. ALT catalyses the transfer of the amino group of L-alanine to α-ketoglutarate to give L-glutamate. The highest levels are found in the liver and the kidneys, and in smaller amounts in heart and skeletal muscle. ALT concentration is increased when hepatic cells are damaged (liver cell necrosis or injury of any cause). Indeed, viral and toxic hepatitis induce a marked elevation of ALT activity in serum. Intake of alcohol, delirium tremens, and administration of various drug induce slight or moderate elevation of ALT. ALT concentration in serum is also slightly increased in various conditions such as: muscular dystrophy, hemolytic disease, myocardial infarction...

ALT is more liver specific than AST (Aspartate aminotransferase). Measurement of both AST and ALT has some value in distinguishing hepatitis from other parenchymal lesions.



ALT serum level can be decreased in case of vitamin B6 deficiency.

METHOD (4)

IFCC method without pyridoxal phosphate (P-5'-P). Kinetic. UV.

PRINCIPLE (4)

Kinetic determination of alanine aminotransferase (ALT) activity:

LDH = Lactate dehydrogenase.

REAGENTS COMPOSITION

Reagent 1 : R1

Tris buffer, pH 7.50 (30°C) 125 mmol/L
L-Alanine 680 mmol/L

LDH ≥ 2 000 U/L
Sodium azide < 0,1 %

Reagent 2 : R2

α-Ketoglutarate 97 mmol/L
NADH 1,1 mmol/L
Sodium azide < 0,1 %

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2 4 x 3 mL
- CONT-0060 ELITROL I 10 x 5 mL
- CONT-0160 ELITROL II 10 x 5 mL

- General Laboratory equipment.

- Normal saline solution (NaCl 9 g/L)

- Do not use materials that are not required as indicated above.

WARNINGS AND PRECAUTIONS

- This reagent kit is for professional *in vitro* diagnostic use only.

- These reagents contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment to avoid contamination.

- Do not interchange reagent vials from different kits.

- For more information, Safety Data Sheet (SDS) is available on request for professional user.

STABILITY OF REAGENTS

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the bottle labels.

On board stability:

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagents are ready to use.

ALT/GPT 4+1 SL

REAGENT DETERIORATION

- The reagent solution should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

DAMAGED PACKAGING

Do not use the reagent if the damages of packaging might have an effect on the product performance (leak-ages, pierced bottle).

SAMPLES (1,5)

Specimen
- Serum and lithium heparinized plasma, free from hemolysis.
- Do not use other specimens.

Warnings and precautions

- According to Good Laboratory Practice, venipuncture should be performed prior to the administration of drugs.

Storage and stability

Samples are stable 3 days at room temperature and 7 days at 2-8 °C. ALT stability is better maintained at -70 °C.

REFERENCE VALUES (1,4)

Men : ≤ 45 U/L
Women : ≤ 34 U/L

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

For ELITech Clinical Systems Selectra Analyzers, applications are available on request.

Wavelength 340 nm
Temperature 37 °C
Read against reagent blank.

Reagent R1	240 µL
Sample	30 µL

Mix and wait 4 minutes and 43 seconds, then add:

Reagent R2	60 µL
------------	-------

Mix and after a 50 second incubation, measure the change of absorbance per minute (ΔA/min.) during 159 seconds.

- With Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

- High ALT values may induce falsely low results due to the depletion of the substrate (total consumption of NADH before reading of the result). For ELITech Clinical Systems Selectra Pro Series Analysers, the application contains a specific alarm to warn the users.

CALCULATION

$$\frac{\Delta A}{\Delta A} \text{ Sample} \times n \quad n = \text{calibrator concentration}$$

$$\frac{\Delta A}{\Delta A} \text{ Calibrator}$$

$$\text{Conversion factor :} \quad U/L \times 0.0167 = \mu\text{kat/L}$$

CALIBRATION

For calibration, multiparametric calibrator ELICAL 2 must be used. Its value is traceable to IFCC reference method⁽⁵⁾.

Calibration frequency : The calibration is specific for each analyser. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

To ensure adequate quality, control sera such as ELITROL I (normal control) and ELITROL II (abnormal control) should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local and legal requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

Measuring range

Determined according to CLSI⁽⁷⁾ EP6-A protocol, the measuring range is from 15.0 to 250.0 U/L (0.25 to 4.17 µkat/L). Samples exceeding 250.0 U/L should be diluted 1:10 in NaCl 9 g/L solution (normal saline) and re-assayed. Use of this procedure extends the measuring range to 250.0 to 2500.0 U/L (4.17 to 41.67 µkat/L).

For Selectra TouchPro users, the « rerun dilution » function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LoD) and Limit of Quantification (LoQ)

Determined according to CLSI⁽⁸⁾ EP17-A protocol.
LoD = 1.8 U/L (0.03 µkat/L)
LoQ = 5.0 U/L (0.08 µkat/L)

Precision

Determined according to CLSI⁽⁹⁾ EP5-A2 protocol.

		Mean	Within-run	Total
	n	U/L	µkat/L	CV (%)
Low level	80	24,6	0,41	2,1
Medium level	80	50,1	0,84	1,6
High level	80	212,4	3,54	0,8

Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-Approved system equipment (IFCC method) on 99 human serum samples according to CLSI⁽¹⁰⁾ EP9-A2 protocol.

The sample concentrations were between 17.3 and 256.4 U/L (0.29 and 4.28 µkat/L).

The parameters of the linear regressions are as follows:
Correlation coefficient: (r) = 0.999
Linear regression: y = 0.963 x - 1.7 U/L (0.03 µkat/L)

Limitations/Interferences

- Hemolysed samples should not be used since significant hemolysis may increase ALT concentration because of high levels of ALT in erythrocytes.

- Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according to EP7-A2 protocol of CLSI⁽¹¹⁾. Recovery within ± 10% of initial value of ALT activity of 40.0 U/L and 200.0 U/L.
Unconjugated Bilirubin : No significant interference up to 30.0 mg/dL (513 µmol/L).
Conjugated Bilirubin : No significant interference up to 29.5 mg/dL (504 µmol/L).
Triglycerides: No significant interference up to 2820 mg/dL (31.86 mmol/L).
Pyruvate: No significant interference up to 3.0 mg/dL.
Ascorbic acid: No significant interference up to 20.0 mg/dL (1136 µmol/L).
Acetylsalicylic acid: No significant interference up to 200.0 mg/dL (11.1 mmol/L).
Acetaminophen: No significant interference up to 30.0 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽¹²⁾

- Many other substances and drugs may interfere. Some of them are listed in Young.⁽¹³⁻¹⁴⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

On board stability/Calibration frequency

On Board Stability: 28 days
Calibration frequency: 28 days
Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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