TRIGLYCERIDES MONO SL NEW

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☛BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Naito, H.K., *Coronary Artery Disease and Disorders of Lipid Metabolism. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C. (Mosby, Inc. eds. St Louis USA), (2003), 603.
- Expert Panel on Detection, *Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP)*. *JAMA*, (2001), **285**, 2486.
- Fossati, P., Prencipe, L., *Serum triglycerides determined colorimetrically with an enzyme that produces hydrogen peroxide. Clin. Chem.*, (1982), **28**, 2077.
- Tietz, N.W., *Clinical guide to laboratory tests*, 3rd Ed. (W.B. Saunders eds. Philadelphia USA), (1995), 610.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), **23** (16).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd edition, AACCC Press (1997).
- Young D.S., *Effects of drugs on clinical laboratory tests*, 4th edition, AACCC Press (1995).

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☛ MATERIALS REQUIRED BUT NOT PROVIDED

- CONT-0060 ELITROL I 10 x 5 mL
- CONT-0160 ELITROL II 10 x 5 mL
- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

☛ WARNINGS AND PRECAUTIONS

- This reagent and this standard are for professional *in vitro* diagnostic use only.
- The reagent R and standard contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these products always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- The standard should be immediately and tightly capped to prevent contamination and evaporation.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

STABILITY OF REAGENTS

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.

On board stability :
 The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagent and standard are ready to use.

☛ REAGENT DETERIORATION

- The reagent and standard solution should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

☛ DAMAGED PACKAGING

Do not use the reagent if the damages of packaging might have an effect on the product performance (leakages, pierced vial).

☛ SAMPLES (4)

- Specimen**
- Serum or lithium heparinized plasma from fasting patients (≥ 12 hours).
- Do not use icteric or hemolyzed samples.
- Do not use other specimens.
- Warnings and precautions**
- According to Good Laboratory Practice, venipuncture should be performed prior to the administration of drugs. A venipuncture could lead to false results if performed during or immediately after the administration of some drugs.
- Collect the samples in tubes and stoppers free of glycerol.
- Separate from cells within 2 hours.
- Storage and stability**
- Samples are stable 5 to 7 days if stored at 2-8 °C, 3 months at -15 to -20 °C and several years at -70 °C. Avoid repeated freezing and thawing.

REFERENCE VALUES (2)

The NCEP (American National Cholesterol Education Program) has established the following classification for triglycerides levels according to the risk of developing coronary heart disease:

Risk Classification Level (mg/dL)	Level (mmol/L)
Normal	< 150 1,69
Borderline high	150 - 199 1,69 - 2,25
High	200 - 499 2,26 - 5,64
Very High	≥ 500 5,65

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

☛ PROCEDURE

For ELITech Clinical Systems Selectra Analyzers applications are available on request
 Wavelength 505 nm
 Temperature: 37 °C
 Read against reagent blank.

	BLANK	CALIBRATION	TEST
Reagent R	300 µL	300 µL	300 µL
Distilled water	3 µL	-	-
Standard	-	3 µL	-
Sample	-	-	3 µL

Mix and read the absorbances (A) after an incubation of 11 minutes and 30 seconds.

- Triglycerides Mono SL New reagent can be contaminated by Cholesterol HDL SL 2G reagent. In order to avoid contamination on Selectra ProM, ProXL, E and XL, program incompatibilities as follows:

Software	Menu	Parameter
TouchPro	Probe	Link / Cholesterol
	incompatibilities	HDL SL 2G -Acid Solution
Other	Needle	Cholesterol HDL SL 2G <<HCl incompatibility

For other instruments, repeat any absurd results after programming a needle wash.

- Lipase SL reagent is strongly contaminated by Triglycerides SL reagent. In order to avoid cuvette contamination on Selectra instruments, program the following incompatibilities:

Software	Menu	Parameter
TouchPro	Test	Link / Triglycerides SL
	incompatibilities	- Acid Solution
Other	Cuvette	Triglycerides SL << HCl incompatibility

In order to avoid needle contamination on Selectra instruments, do not program Lipase SL and Triglycerides SL in the same run. Ensure the instrument goes back to "stand-by" status before launching a run containing Lipase SL.

CALCULATION

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times n \quad n = \text{standard concentration}$$

Conversion factor : mg/dL x 0.0113 = mmol/L
 mg/dL x 0.01 = g/L

CALIBRATION

Concentration value for Triglycerides Standard 200 mg/dL is traceable to the reference method ID-MS (Isotope Dilution - Mass Spectrometry).

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

☛ QUALITY CONTROL

To ensure adequate quality, control sera such as ELITROL I (normal control) and ELITROL II (abnormal control) should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local and legal requirements.

☛ PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

- **Measuring range**
 Determined according CLSI EP6-A protocol⁽³⁾, the measuring range is from 30 to 1000 mg/dL (0.34 to 11.30 mmol/L).
 Samples exceeding 1000 mg/dL should be diluted 1:5 with NaCl 9 g/L solution (normal saline) and re-assayed. Use of this procedure extends the measuring range to 1000 to 5000 mg/dL (11.30 to 56.50 mmol/L).

For users with Selectra TouchPro software, the «rerun dilution» function performs the sample dilution automatically. Results take the dilution into account.

- Precision

Determined according to CLSI EP5-A2 protocol⁽³⁾.

		Mean	Within-run	Total
	n	mg/dL	mmol/L	CV (%)
Low level	80	44	0.50	2.0 3.8
Medium level	80	131	1.48	0.9 2.3
High level	80	267	3.02	1.2 2.4

- Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (enzymatic and colorimetric method) on 99 human serum samples according to CLSI EP9-A2 protocol⁽⁷⁾.

The sample concentrations were between 30 and 957 mg/dL (0.34 and 10.81 mmol/L). The parameters of the linear regressions are as follows :
 Correlation coefficient: (r) =0.999
 Linear regression: y = 1.019 x + 1 mg/dL (0.10 mmol/L)

- Limitations/Interferences

- Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according CLSI EP7-A2 protocol⁽³⁾. Recovery is within ±10% of initial value of Triglycerides concentration of 133 mg/dL and 266 mg/dL.

Unconjugated Bilirubin: No significant interference up to 15 mg/dL (257 µmol/L).

Conjugated Bilirubin: No significant interference up to 5.9 mg/dL (101 µmol/L).

Hemoglobin : No significant interference up to 125 mg/dL.

Uric acid : No significant interference up to 24.2 mg/dL (1440 µmol/L).

Ascorbic acid : No significant interference up to 2.0 mg/dL. Concentrations above the therapeutic levels will interfere and cause erroneous results.

Methyl-dopa : No significant interference up to 1.0 mg/dL (47 µmol/L).

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽⁹⁾

- Results can be falsely lowered by significant levels in sample of NAC (N-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metamizole.

- Many other substances and drugs may interfere. Some of them are listed in Young.⁽¹⁰⁻¹¹⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

- On board stability/Calibration frequency

On Board Stability: 28 days
 Calibration frequency: 14 days
 Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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